



Ask for quote

## COVID-19 Rapid Test (Nasal Swab)

Covid-19 test

Rapid test  
Covid-19

**Newgene Covid-19 Rapid Test Nasal (Swab)** approved as self-test in some countries. It is an in vitro diagnostic test for the qualitative detection of novel coronavirus antigens in Nasal swab, using the rapid immunochromatographic method. The identification is based on the monoclonal antibodies specific for the novel corona virus antigen. It will provide information for clinical doctors to prescribe correct medications. The advantages speak for themselves: The packages are available in 3 different sizes: 1-pack, 5-pack & 25-pack

- 10 minutes for results
- 99,2% Specificity
- 97,1% Sensitivity
- Easy to operate
- Test can be carried out everywhere

[Download Product Information](#)

**Varumärke**

**Leveranstid**

2-3 work weeks

**Forpackning**

1-pack

25-pack

5-pack