



Ask for quote

## Disposable surgical mask Type IIR sterile with tie bands

Facemask

Sterile liquid repellent surgical mask. Masks in three layers for single use in healthcare. Adjustable nose clip and approx 37 cm tie bands. CE marked and approved according to EN14683: 2019. ISO 22609:2004. 20 pcs / box.

### Farg

Blue

### Material

Nonwoven PP  
+ Filter PP +  
Nonwoven PP

### Produktens-hojd

95 mm

### Produktens-breidd

175 mm

### Produktens-vikt

3 g

### Box

Yes

### Polybag

No

### Statistik-nummer

6307909100

### Antal-ytterkartong

2500 pcs

### Kartongens-storlek

65 x 47 x 39

cm

### Kartong-bruttovikt

12 kg

### Leveranstid

Approx 3 work

weeks

### Ursprungsland

China