



Ask for quote

Filtering half mask FFP3 Cup

Facemask
FFP3, filtering half mask to protect against particles. Disposable cup shaped protective face mask. Tested and approved by notified body according to EN149:2001+A1:2009.
Farg
White
Material
4 layers (Acupuncture cotton, Melt-blown fabric, PTFE fabric, Acupuncture cotton)
Certifikat
Produktens-hojd
12,7 (+/- 0,5) cm
Produktens-bredd
15,2 (+/- 0,5) cm
Box
Yes
Antal-ytterkartong
360 st
Kartongens-storlek
47,5 x 42 x 49,5 cm
Kartong-bruttovikt
9,1 kg
Leveranstid
Approx 4 work weeks
Statistik-nummer
6307909100
Ursprungsland
China