



Ask for quote

COVID-19 Rapid Test Saliva

Covid-19 test

Rapid test
Covid-19

Newgene Covid-19 Rapid test (saliva) is an in vitro diagnostic test for the qualitative detection of novel coronavirus antigens, using the rapid immunochromatographic method. The identification is based on the monoclonal antibodies specific for the novel corona virus antigen. It will provide information for clinical doctors to prescribe correct medications. The advantages speak for themselves: The packages are available in 3 different sizes: 1-test, 5-tests & 25-tests

- 10 minutes for results
- 99,2% Specificity
- 97,1% Sensitivity
- Easy to operate
- Test can be carried out everywhere

[Download Product Information](#)

Varumärke

Leveranstid

2-3 work weeks

Forpackning

1-pack

25-pack

5-pack