



Ask for quote

Retail Surgical mask Type IIR

Retail

SkyScreen® surgical mask Type IIR, a three ply disposable mask with high bacterial filter efficiency. A face mask can protect the surroundings from microorganisms in your exhaled air and protect you from drops and splashes. Has a low respiratory resistance to be easy to breathe through. SkyScreen® surgical mask Type IIR consists of three filtration layers, polypropylene fiber fabric, meltblown filter, nose clip and elastic ear straps. Does not contain latex. EN 14683:2019 + AC:2019 The package contains 10 face masks. EU representativeMEC Medicare

Color



Material

Nonwoven PP
+ Filter PP +
Nonwoven PP

Produktens-hojd

95 mm

Produktens-vikt

3 g

Box

Yes, for retail

Polybag

No

Statistik-nummer

6307909100

Antal-innerkartong

10 pcs

Antal-ytterkartong

2000 pcs

Kartong-bruttovikt

10 kg

Minsta-antal

25000 pcs

Leveranstid

Approx 3 work
weeks

Ursprungsland

China